



## Why FloorBed™?

The lowest bed in its class currently available, FloorBed™ is just 7.1cm high - similar to a tennis ball. With a low profile mattress giving a combined height of only 15cm it significantly reduces the risk of injury if users roll out of bed. FloorBed™ rises to full nursing height for care duties.

## Functions & Benefits

### Benefits of the low rise function

- Reduces risk of injury following a fall from bed
- Reduces requirement for bed rails
- Removes risks associated with bed rails
- Reduce user agitation
- Reduces requirement for crash mats
- Independent transfers
- Independent leg lift into bed

### Benefits of profiling bed functions

- Full Nursing height
- Compatible with specialist mattress systems
- No trailing wires
- Back rest/elevation
- Knee break
- Pressure area care position
- Bed Extension

“The fact that the bed lowers to the floor ensures that a resident cannot fall out, but can roll out of the bed. When they are at floor they are usually unable to rise unaided from the floor so will not wander or fall in their room.”

HOME MANAGER

“We put the Resident in the FloorBed™ and it was really beneficial. The risk of injury to her was immediately reduced and if she rolled out there was no longer any distance for her to fall. This also reduced the time staff had to monitor her as the risk of injury was so greatly reduced. During the time we had the bed the lady rolled out a number of times but sustained no injuries whatsoever.”

HOME MANAGER

“Excellent bed for mobile residents with dementia or reduced mobility.”

HOME MANAGER



[www.floorbed.com](http://www.floorbed.com)



### Manufactured by:

Accora Ltd  
 Charter House  
 Barrington Road  
 Orwell  
 Cambridge, SG8 5QP  
 T: 01223 206100  
 F: 01223 206120  
 info@accora.uk.com  
 www.accora.uk.com  
 Registered to ISO 9001:2008



### Reduce Injury

Total bed and mattress height of FloorBed™ is just 15cm - a bed fall becomes a roll, almost eliminating the risk of injury

### Reduce Carer Burden

Reduced workload and strain for carers as less supervision is needed

### Manage Risk

Safeguard your organisation by mitigating exposure to user and corporate risk



Total height of 15cm when used with LowMatt low profile mattress



No trailing wires



Full nursing height




Concealed wheels



Backrest and kneebreak functions



Foldable for easy transport



Mattress platform height of just 7.1cm – similar to a tennis ball!

#### TECHNICAL SPECIFICATION

Dimensions:	92 x 224cm (3ft x 7ft 4in)	Lifting range: (excl. mattress)	7.1cm - 64cm (3in - 2ft)
Platform size:	90 x 200cm (3ft x 6ft 7in)	Profiling functions:	Backrest, Kneebreak, Anti-Trendelenburg, Trendelenburg
Weight Capacity:	180kg (28st)	Part number:	NSB-0-FL1-100
Max user weight: (with mattress weight of 10kg)	170kg (27st)		